# **ThermoML: User Guide**

ThermoML is a thermal image-based system for detecting inflammation in hands using AI models and computer vision. It features a web interface, Firebase-based backend, and a Python-powered inference system with support for local and Hugging Face Space deployment.

Live Demo: <a href="https://thermoml.web.app">https://thermoml.web.app</a> Contact: <a href="mailto:ThermoML@outlook.com">ThermoML@outlook.com</a>

> Note: Backend image processing is currently offline. To request temporary access to the server, please contact us via email.

# **Getting Started:**

## **Clone the repository:**

git clone https://github.com/ShaharZeharia/TheremoML.git cd TheremoML

## **Backend**

### requirements:

- Python 3.12.0
- Git
- ngrok(optical, for public url)

## **Download required models:**

- <u>SAM</u> facebookresearch/segment-anything Checkpoint: vit h
- Google Download Latest checkpoint
- Our built Resnet50 models Google Drive Link

#### **Create .env for the Python App:**

GOOGLE\_APPLICATION\_CREDENTIALS=path/to/your/credentials.json MODELS\_PATH=./models/ PORT=8000 DEBUG=True CALLBACK\_SECRET=your\_secure\_callback\_secret

## Run the Python Backend:

cd app pip install -r requirements.txt python main.py (Optinal) Expose with ngrok: ngrok http 8000

# Firebase Cloud Functions(Local)

**Create .env inside /functions** 

HUGGINGFACE\_URL=https://your-ngrok-url.ngrok.app/api HUGGINGFACE\_TOKEN=your\_huggingface\_token (or leave blank if local) CALLBACK\_KEY=your\_secure\_callback\_secret #Same as the one in the python .env file

#### **Install and Run:**

npm install -g firebase-tools firebase login firebase use --add cd functions npm install firebase emulators: start --only functions

## **Frontend**

The frontend React app is located in the root folder. Install Frontend Dependencies: npm install
Run Frontend Locally: npm run
This will launch the app at http://localhost: 3000.

# **Hosting**

Build the react app:

npm run build

Then deploy the build/folder.

## **Contact**

If you encounter issues or want to contribute, feel free to reach out:

ThermoML@outlook.com